

# Introducing the *Flow-Safe*<sup>™</sup> CPAP System

The First Disposable CPAP System  
with built-in Safety Features for less.



Disposable, Latex-Free  
*Flow-Safe* CPAP System

- Built-in Manometer
- Built-in Pressure Relief Valve
- Less Oxygen Consumption

For the treatment  
of Acute Pulmonary Edema (APE).

# The first disposable CPAP System for Acute Pulmonary Edema with built-in manometer and pressure relief valve

## Flow-Safe™ CPAP Systems



## Flow-Safe™

• <i>Built-In Manometer on Every Device</i> Manometer Accuracy +/- 3 H <sub>2</sub> O up to 15 cm H <sub>2</sub> O.	X
• <i>Built-In Pressure Relief Valve</i> Limits Maximum CPAP Pressure to 25 cm H <sub>2</sub> O +/- 5 cm H <sub>2</sub> O @ 25 lpm.	X
• <i>Highly Portable (No costly flow generator repairs or downtime)</i> Completely Disposable, Latex-Free CPAP System.	X
• <i>Easy Set-Up</i> Includes Mask, Manometer, Pressure Relief Valve and Head Harness, just requires O <sub>2</sub> source and flowmeter to operate. CPAP unit only also available.	X
• <i>Color-Coded for Easy Size Identification.</i>	X
• <i>Less Oxygen Consumption</i>	X
• <i>Allows Suctioning</i> Without Removing the Mask.	X
• <i>Clinicians can add a Nebulizer</i> Combining Aerosol Delivery and CPAP Therapies.	X
• Accepts 22mm Bacterial Filter.	X
• Accepts 15 mm/22 mm CO <sub>2</sub> Sampling Lines without additional adapters.	X
• Standard 15 mm/22 mm port size.	X
• Choice of mask styles	X

# Flow-Safe™

Soft Cushioned Mask with built-in manometer and pressure relief valve



Deluxe Contoured Mask with built-in manometer and pressure relief valve

## Product Specifications

Lightweight and Portable	75 grams nominal (less mask and harness) 61 mm X 60 mm X 50 mm (unit only)
<u>Flow (LPM)</u>	<u>CPAP/PEEP (cm H<sub>2</sub>O)</u>
10	1.5 - 2.0
15	3.0 - 4.0
20	6.0 - 7.0
25	8.5 - 10

## Ordering Information

## Quantity

#10-57000	Flow-Safe CPAP System with large adult soft cushion mask, manometer, pressure relief valve, head harness, oxygen tubing (7')	5/Box
#10-57001	Flow-Safe CPAP System with small adult soft cushion mask, manometer, pressure relief valve, head harness, oxygen tubing (7')	5/Box
#10-57002	Flow-Safe CPAP System with child soft cushion mask, manometer, pressure relief valve, head harness, oxygen tubing (7')	5/Box
#10-57003	Flow-Safe CPAP System Only with manometer and pressure relief valve, oxygen tubing (7')	5/Box
#10-57006	Flow-Safe CPAP System with large adult deluxe contoured mask, manometer, pressure relief valve, head harness, oxygen tubing (7')	5/Box
#10-57007	Flow-Safe CPAP System with small adult deluxe contoured mask, manometer, pressure relief valve, head harness, oxygen tubing (7')	5/Box
#10-57008	Flow-Safe CPAP System with child deluxe contoured mask, manometer, pressure relief valve, head harness, oxygen tubing (7')	5/Box
#10-57100	Large adult deluxe contoured mask	5/Box
#10-57101	Small adult deluxe contoured mask	5/Box
#10-57102	Child deluxe contoured mask	5/Box
#10-57004	Disposable Head Harness	20/Box
#10-57005	Silicone Head Harness	1/Ea.
#10-55376	Bacterial Filter (22 mm)	20/Box

The first disposable CPAP System with built-in manometer and pressure relief valve.

*Highly Portable, Easy Set-Up, Built-In Manometer & Pressure Relief Valve*